
Patient Instruction: Urinary Stone Collection (CalCuCatch)

What are patient instructions?

Patient instructions are written toward helping explain the proper way in collecting laboratory specimens. If followed correctly, you are ensuring your test results will be accurate and reliable. They are instructions to be followed before, during and after the collection process. This is a vital step in the outcome of your laboratory results.

Contacts:

If you have any questions before, during or after the collection, feel free to contact the following:

Nekoosa Clinic Laboratory
Monday-Friday; 8 a.m. to 5 p.m.
(715) 886-2140

Lakes Clinic Laboratory
Monday-Friday; 8:30 a.m. to 5 p.m.
(715) 325-8340

Wisconsin Rapids Clinic Lab Area
Monday-Friday; 8 a.m. to 5 p.m.
(715) 422-7767

Riverview Hospital - Main Laboratory
Monday-Friday; 7 a.m. to 7 p.m.
Saturday-Sunday; 8 a.m. to 3 p.m.
(715) 421-7433

Instructions:

Your doctor has ordered a urinary stone analysis.

- 1) Be sure to pick up a CalCuCatch “Urinary Stone Interceptor” from the clinic or lab.
- 2) Remove the lid and set it aside. We will need this in Step 6.
- 3) Open the collection cup to its fullest position (telescope) and begin urinating directly into the cup. Notice that a screen is attached to the bottom. Be sure you are over a toilet when performing this step.
- 4) If you look at the lid, you will notice that there is a compartment in the middle. Open the smaller compartment lid and place any material that collected on the screen after urinating into the compartment. Replace the smaller lid.
- 5) Rinse the cup and dry. Break down the cup back to its original position. Replace the larger lid.
- 6) Place the interceptor back into the ziplock bag. Drop off the specimen at the clinic or lab location as specified.

Specimen Delivery:

Keep the specimen at room temperature after collection and drop off at the clinic or lab location as soon as possible as specified. Refer to the contact information for hours of operation.

Handwashing:

It is very important to ensure proper handwashing is achieved after the collection process. Wash hands for 3-5 minutes with soap and water, rinse thoroughly.

Identification:

Be sure the CalCuCatch is labeled with your name, date of birth, date of collection, and time of collection.

Our Commitment

We are committed to help and assist you in every possible way. This is done by providing the most accurate and precise patient instructions. We hope your experience will be positive.

If you feel we can improve these instructions, please call or email the Clinic Laboratory Supervisor directly at any time. We value your input in making our services better.

Riverview Family Clinic
Laboratory Supervisor
(715) 886-2141

Email: hawmic@rhahealthcare.org

Lab Appointment Needed:

Please call ahead and schedule your lab appointment if you need lab work done or dropping off your specimen.

Nekoosa Clinic
(715) 886-2100

Lakes Clinic
(715) 325-8300

Wisconsin Rapids Clinic
(715) 421-7474

Thank You

Riverview Family Clinic

Nekoosa Clinic

Nekoosa Medical Center
1015 Angelus Drive
Nekoosa, WI 54457
715-886-2100

Lakes Clinic

Lakes Area Medical Center
1160 Rome Center Drive
Nekoosa, WI 54457
715-325-8300

Wisconsin Rapids Clinic

Riverview Medical Center
410 Dewey Street
Wisconsin Rapids, WI 54494
715-421-7474

Web Site:

www.riverviewclinic.net

Urinary Stone Collection (CalCuCatch)

Clinical Laboratory Patient Instructions

